




PAD Below-The-Knee: Enhancing the Performance of Peripheral Vascular Catheters

Advancements in minimally invasive therapies have revolutionized the treatment of peripheral artery disease below-the-knee (BTK). However, despite advances, device makers face critical challenges in developing newer, more capable devices to meet growing market demand.



Market: Medical Device

Sub-Market: Peripheral

Body Part: Leg

Condition: Peripheral Artery Disease (PAD)

Device: Peripheral Vascular Catheters

Zeus Product: StreamLiner™ NG



Peripheral Artery Disease



Peripheral artery disease (PAD) is the narrowing of the blood vessels in the arms or legs caused by atherosclerosis – a buildup of plaque in and on the walls of an artery. Left untreated, PAD can progress over time, leading to serious complications, including heart attack, stroke, and critical limb ischemia (CLI).

PAD is a growing global health concern. While figures vary, PAD is estimated to affect between 200–236 million people worldwide.^[1,2] When looking at trends in the prevalence of PAD, systematic reviews have estimated the number of people living with PAD increased by roughly 45% between 2000 and 2015.^[2,3,4] It's clear that despite advancements in medicine around the world, PAD remains a growing health concern.

The treatment of PAD largely depends on the severity of the disease. It can range from simple lifestyle changes to open surgery or, in the most severe cases, amputation of the affected limb. However, advancements in endovascular devices and surgical techniques have resulted in more widespread adoption of minimally

invasive surgical interventions to treat this condition.

From guide to delivery catheters, there are a multitude of devices used in these minimally invasive interventions, especially in treating lower extremity PAD – the most common type. As the prevalence of PAD has continued to rise, so too has demand for catheters capable of delivering new, innovative therapies into the complex vasculature below-the-knee (BTK), with industry reports suggesting the global peripheral vascular devices market is projected to grow at a compound annual growth rate (CAGR) of 5.8% from 2025 to 2030.^[5]

200 Million

People Affected By PAD

Source: <https://www.ahajournals.org/doi/10.1161/CIR.0000000000001153>

Global peripheral vascular devices market projected to expand at a CAGR of

5.8%

Source: <https://www.grandviewresearch.com/industry-analysis/peripheral-vascular-devices-market>



Below-The-Knee Constraints

While many options exist for treating peripheral artery disease, ranging from lifestyle changes and medication to open surgery, the most common approach is via minimally invasive catheter-based procedures.

In these procedures, a series of catheters are guided through the vasculature to the location of the arterial blockage with the ultimate goal of restoring blood flow.

Perhaps the most critical aspect of these procedures is successfully navigating the catheters to and from the treatment site without failure. Regardless of the interventional approach, the tortuous vasculature below-the-knee requires peripheral vascular catheters to not only be incredibly soft and flexible, but also possess high tensile and burst strength in order to successfully deliver therapies without failure – a difficult engineering feat.

Furthermore, as the need for even smaller, longer, and more capable catheters continues to rise, device engineers are faced with progressively more difficult material selection and design dilemmas, often becoming increasingly reliant on the catheter liner itself to shoulder much of the burden of the ever-evolving performance requirements of next-gen designs.

While many catheters used in the treatment of PAD utilize film-cast PTFE liners due to their desirable mix of ultra-thin walls, lubricity, and flexibility – there are limitations. Inherently, existing film-cast liners are prone to surface imperfections and pinholes as a byproduct of the legacy casting process, which can have a negative impact on mechanical performance and reliability, as well as manufacturing yield. As such, for continued device innovation in this field, it is clear that more advanced film-cast PTFE liners are needed.



Top Design Challenges & Considerations

- Small Diameter Vasculature
- Tortuous and Delicate Anatomical Pathways
- Need for Soft and Flexible, Yet Strong Catheter Liners



Minimize Defects, Maximize Strength and Performance
 INTRODUCING STREAMLINER™ NG

An ultra-thin and flexible film-cast PTFE liner engineered for *maximum strength* with *minimal defects*.

As a result of Zeus' proprietary film-cast process, which minimizes defects and pinholes, StreamLiner™ NG catheter liners take flexibility, mechanical performance, and reliability to the next level while still featuring the exceptional sizing and tolerances that the StreamLiner™ series is known for.

Remarkable Burst & Yield Strength

Zeus' proprietary film-cast process results in remarkable burst and yield strength, helping enhance mechanical performance of the finished device.

Minimal Scrap

Fewer pinholes, fewer problems. The ultra-low occurrence of imperfections and pinholes in StreamLiner™ NG liners can help improve manufacturing efficiencies and production yields.

High Concentricity

StreamLiner™ NG liners feature highly consistent wall thicknesses, helping to provide tighter OD tolerances for the finished device, as well as improved consistency in device performance.

Exceptional Sizing & Flexibility

True to the StreamLiner™ name, StreamLiner™ NG catheter liners are available in a wide range of sizes, with IDs as low 0.017" (0.432 mm) and nominal wall thicknesses down to 0.0005" (0.0127 mm).



FEATURES

StreamLiner™ NG

Proprietary Film-Cast Process

Zeus' novel film-cast process produces exceptionally consistent liners with minimal defects and pinholes, helping pave the way for lower manufacturing costs and higher production yields.

Remarkable Strength

The remarkable burst and yield strength of StreamLiner™ NG provides engineers with more options for designing devices capable of successfully delivering therapies below-the-knee.

Extremely Thin & Uniform Walls

Nominal wall thicknesses as low as 0.0005" (0.0127 mm). These ultra-thin walls enable catheter designs to maintain a minimal outer diameter to successfully traverse the complex vasculature of the lower extremities, or a maximized inner diameter for improved delivery performance of therapies below-the-knee.

Tight Tolerances

Wall thickness tolerance of $\pm 0.00025''$ (0.00635 mm). Our state-of-the-art processes allow us to manufacture film-cast PTFE liners to tight tolerances, enabling device engineers to create catheters to exact specifications for use in specialized applications.

Biocompatible

StreamLiner™ liners are manufactured from USP Class VI materials, giving OEMs important assurance that these liners are safe for use within the body.

High Lubricity

PTFE's low coefficient of friction allows devices, such as drug-coated balloons (DCBs), to slide through the delivery system more easily, helping them reach the treatment location without failure.





AVAILABLE SIZES


StreamLiner™ Series

All PTFE StreamLiner™ catheter liners are produced based on customer specifications, and the table below is a general capability guide. All material can be etched on the OD to enhance adhesion. For even greater jacket-to-liner adhesion, a tie layer coating may also be applied.

PTFE StreamLiner™			
	StreamLiner™	StreamLiner™ OTW	StreamLiner™ NG
MATERIAL	PTFE	PTFE	PTFE
PROCESS	Free-Extruded	Extruded Over-The-Wire	Proprietary Film-Cast
MANDREL	None	Silver-plated copper, Stainless steel	Silver-plated copper, Stainless steel
INSIDE DIAMETER (ID)	0.004" - 0.120" (0.102 mm - 3.048 mm)	0.013" - 0.0915" (0.330 mm - 2.3241 mm)	0.017" - 0.0915" (0.432 mm - 2.3241 mm)
ID TOLERANCE	± 0.0005" - 0.001" (± 0.0127 mm - 0.025 mm)	± 0.0005" (± 0.0127 mm)	± 0.0005" (± 0.0127 mm)
NOMINAL WALL THICKNESS	0.0005" - 0.00075" (0.0127 mm - 0.01905 mm)	0.0004" - 0.00075" (0.0102 mm - 0.01905 mm)	0.0005" - 0.00075" (0.0127 mm - 0.01905 mm)
WALL TOLERANCE	± 0.00025" (± 0.00635 mm)	± 0.0002" - 0.00025" (± 0.0051 mm - 0.00635 mm)	± 0.00025" (± 0.00635 mm)
CUT LENGTH	86" Max.* (2184.4 mm Max.)	86" Max. (2184.4 mm Max.)	86" Max. (2184.4 mm Max.)
SURFACE TREATMENTS	Etched, Tie Layer	Etched	Etched, Tie Layer
STERILIZATION METHODS	Autoclave, EtO	Autoclave, EtO	Autoclave, EtO
STRENGTH	●●●●○	●●●●○	●●●○
FLEXIBILITY	●●●●○	●●●●○	●●●●●

*Liners with a Tie Layer have a maximum cut length of 78" (1981.2 mm).

StreamLiner™ and StreamLiner™ OTW may be shipped with product labels that feature VT, XT, or UT size designations. VT represents a standard nominal wall thickness of 0.00075" (0.01905 mm). XT represents a standard nominal wall thickness of 0.0005" (0.0127 mm). UT represents a standard nominal wall thickness of 0.0004" (0.0102 mm)



Advancing PAD Outcomes

Lower extremity PAD is one of the most under-diagnosed atherosclerotic diseases, resulting in severe pain and limited physical mobility for patients, as well as a significant risk of amputation and associated increase in mortality.

As the prevalence of PAD continues to grow worldwide, so too will the demand for minimally invasive procedures to address the potentially fatal condition.

While minimally invasive catheter-based procedures remain an effective treatment option for peripheral artery disease, the desire to improve the efficacy in harder-to-reach diseased areas of the lower extremities remains a critical challenge for medical device OEMs – one that requires advanced catheter componentry to achieve.

Existing film-cast PTFE liners face inherent strength limitations as a byproduct of the legacy casting process, potentially limiting their effectiveness in certain applications. StreamLiner™ NG catheter liners help engineers overcome these limitations.

By embracing StreamLiner™ NG, engineers can design stronger and more effective peripheral vascular catheters capable of successfully delivering these critical limb- and life-saving therapies.



References

[1] Allison MA, Armstrong DG, Goodney PP, et al. Health Disparities in Peripheral Artery Disease: A Scientific Statement From the American Heart Association. *Circulation*. 2023;148(3). Accessed March 7, 2024. doi:<https://doi.org/10.1161/cir.0000000000001153>

[2] Aday AW, Matsushita K. Epidemiology of Peripheral Artery Disease and Polyvascular Disease. *Circulation Research*. 2021;128(12):1818-1832. Accessed March 7, 2024. doi:<https://doi.org/10.1161/circresaha.121.318535>

[3] Song P, Rudan D, Zhu Y, et al. Global, regional, and national prevalence and risk factors for peripheral artery disease in 2015: an updated systematic review and analysis. *The Lancet Global Health*. 2019;7(8):e1020-e1030. Accessed March 7, 2024. doi:[https://doi.org/10.1016/s2214-109x\(19\)30255-4](https://doi.org/10.1016/s2214-109x(19)30255-4)

[4] Fowkes FGR, Rudan D, Rudan I, et al. Comparison of global estimates of prevalence and risk factors for peripheral artery disease in 2000 and 2010: a systematic review and analysis. *The Lancet*. 2013;382(9901):1329-1340. Accessed March 7, 2024. doi:[https://doi.org/10.1016/s0140-6736\(13\)61249-0](https://doi.org/10.1016/s0140-6736(13)61249-0)

[5] Peripheral Vascular Devices Market Size & Share Report, 2030. Grand View Research. Accessed March 7, 2024. <https://www.grandviewresearch.com/industry-analysis/peripheral-vascular-devices-market>



zeusinc.com | info@zeusinc.com

Americas: +1 803 268 9500 | Europe: +353 74 9109700 | Asia-Pacific: +(86)20 38254906